

# IBC MANAGEMENT SOFTWARE

WORKFLOW MANAGEMENT SOFTWARE FOR NUCLEAR MEDICINE & RADIOPHARMACY



USER EXPERIENCE | TRACEABILITY OF DATA | COMPLIANCE



**COMECER**

AN **ATS** COMPANY

# IBC MANAGEMENT SOFTWARE

WORKFLOW MANAGEMENT SOFTWARE FOR NUCLEAR MEDICINE AND RADIOPHARMACY

## Unlock the power of Streamlined Nuclear Medicine and Radiopharmacy operations!

*“Experience seamless workflow management software, meticulously designed to meet the unique needs of your nuclear medicine and radiopharmacy departments. Our comprehensive solution ensures the precise delivery of the right radiopharmaceutical dose to each patient, at the right time.”*

Effortless Workflow Management without paper!

Effortlessly manage your operations with our software. Seamlessly track and trace all processes within your facility, ensuring secure and transparent data management. Say goodbye to manual, paper-based processes and hello to efficiency.

Total data traceability

Achieve full control and traceability over all aspects of your workflow. Our system provides seamless integration with Comecer hardware, allowing you to monitor and trace each step securely. From radiopharmaceutical preparation, to ordering stock, patient injections and waste monitoring, you'll have a comprehensive overview at your fingertips.

Stay compliant with ease

Rest easy knowing our system is fully compliant with regulatory standards. It simplifies compliance processes, reduces the risk of errors or penalties, and ensures adherence to strict regulations.

*Features you can expect. The features set depends on what variant of IBC you choose.*

*Let our sales team help you choose.*

### Workflow Efficiency and Patient Safety

- Easily guide the preparation of (patient) doses using customizable protocols.
- Record precise administration details for each patient, ensuring safety and compliance.
- Seamlessly import worklists from Hospital Information Systems for efficient dose preparation.
- Monitor user activities within the system for enhanced accountability.

### Stock Management and Traceability

- Track the storage and usage of stock items with real-time visibility.
- Streamline order management and maintain accurate stock records.
- Monitor and manage ordered multidose vials and syringes effectively.
- Handle expired stock by allocating it to designated waste containers, ensuring proper disposal.

### Quality Control and Compliance

- Conduct built-in pharmaceutical quality checks, including Mo-99 Breakthrough Tests and more.
- Perform activity measurements with supported dose calibrators for accurate dosing.
- Customize and define your own dose calculation methods for precise calculations.
- Create and execute preparation protocols to maintain consistency and compliance.

### Integration and Injection Tracking

- Seamlessly integrate with the Comecer IRIS system to accurately track injections.
- Design and print labels with ease for proper identification and labeling.
- Monitor the current nuclear activity in your lab for efficient resource management.


Available solutions

IBC Management Software offers **three powerful variants** of our comprehensive workflow management software, each tailored to meet the unique needs of your department.

## IBC CLINIC

### DESIGNED FOR NUCLEAR MEDICINE DEPARTMENTS


IBC Clinic focuses on patient and dose management within nuclear medicine departments. This variant comes equipped with our standard features for efficient workflow management. Additionally, it includes specialized features that enable seamless integration with hospital information systems, optimizing collaboration and streamlining the preparation and injection of patient doses.



## IBC NM

### EMPOWERING NUCLEAR MEDICINE DEPARTMENTS

IBC Nuclear Medicine is specifically designed to meet the stringent demands of nuclear medicine departments. It encompasses all our standard features for comprehensive workflow management, and goes a step further by offering specialized features. These include end-to-end traceability, seamless integration with Comecer hardware, and an Extended Protocol and Quality Control for pharmaceuticals functionality. With IBC Nuclear Medicine, you can standardize your workflow, optimize efficiency, and deliver exceptional patient care.



## IBC RP

### TAILORED FOR RADIOPHARMACIES

IBC Radiopharmacy is a dedicated variant for radiopharmacies, ensuring smooth operations in the production workflow of radiopharmaceuticals. It encompasses our standard features for efficient workflow management, while offering specialized functionalities to meet the specific needs of radiopharmacies. With IBC Radiopharmacy, you can take control of your radiopharmacy operations, ensuring compliance with regulatory standards, streamlining processes, and enhancing overall efficiency.



## Ready to transform your operations?

Contact our sales team today for a free consultation and personalized guidance. We'll determine the best solution for your department and provide a comprehensive software demo. Receive a detailed quotation tailored to your needs and get assistance with installation and validation, ensuring a smooth transition.

**We understand that you may have some concerns, and we're here to address them.**

### Cost

Rest assured that our system offers long-term cost savings, improved efficiency, and reduced errors, providing a positive return on investment.

### Complexity

Our user-friendly software ensures a seamless transition, backed by comprehensive support, training, and onboarding services.

### Compliance

Our system is fully compliant with regulatory standards, streamlining compliance processes and mitigating risks.

### Customization

We offer a flexible solution that can be customized to align with your specific needs and workflows.

### Integration

Our system is designed for seamless integration with common (Comecer and some third-party) equipment and systems used in nuclear medicine departments.

Unlock the power of streamlined operations  
with our comprehensive software solution and deliver exceptional patient care.  
**We hope to speak to you soon!**


WE ARE YOUR PARTNER FOR LIFE

#### HEADQUARTERS

##### COMECER S.p.A.

Via Maestri del Lavoro, 90  
48014 - Castel Bolognese (RA) - Italy  
t: +39 0546 656375 - f: +39 0546 656353  
comecer@comecer.com - www.comecer.com



 Legal manufacturer:  
Comecer Netherlands B.V.  
Madame Curieweg 1  
8501 XC Joure - NL

#### SUBSIDIARY COMPANIES

##### COMECER Inc.

5521 Meadowbrook Industrial Court  
Rolling Meadows - IL 60008 - USA  
t: +1 (847) 925-1238  
serviceamerica@comecer.com

##### COMECER Netherlands B.V.

NL-8501 XC - Joure - The Netherlands  
t: +31 (0) 513 41 69 64 - f: +31 (0) 513 41 69 19  
info\_nl@comecer.com

##### COMECER India Pvt Ltd

Office No.112-113, Veena-Beena Arcade, Gurunanak Road  
Bandra (West) - Mumbai - India - 400050  
t: (+91) 22 26423421 - f: (+91) 22 26423412



www.comecer.com