

A black medical device, the Klarity Athena Supine Breast System, is shown from a top-down perspective. It features a central breast support area with two large, rounded handles on either side. The device is mounted on a long, narrow base with a cushioned seat at the bottom. The Klarity logo is visible on the side of the base. The image is set against a white background with decorative green and pink lines on the left side.

Klarity[®]

ATHENA

**SUPINE BREAST
SYSTEM**

ATHENA™ SUPINE BREAST SYSTEM

The ultimate versatility for precise, stable supine breast and torso treatments

The Klarity Athena™ Supine Breast System offers advanced solutions for precise supine breast and torso positioning. The Athena™ allows for comprehensive arm positioning possibilities with arm cups, upper hand grips, and lower hand grips for asymmetrical positioning. The included Wingboard can be configured with T-bar or U-grip Handles to support customizable shoulder placement. The system also features a securely locking, easily adjustable bum-stop for enhanced stability.

Two Athena™ models are available, each with distinct elevation angles for seamless clinical integration. Both the Athena™ baseplate and Wingboard are compatible with standard S-type head masks, and can also be used with the Klarity SRS Adaptor. Specialized breast and torso masks are also available to meet diverse treatment requirements.

- Extremely lightweight
- Comprehensive arm placement options
- Compatible with S-type head masks, as well as breast and torso masks



R610-PDCF5	R610-PDCF6	R610-PDKF5	R610-PDKF6
Carbon Fiber	Carbon Fiber	Kevlar MR	Kevlar MR
145cm x 41cm	145cm x 41cm	145cm x 41cm	145cm x 41cm
6.3kg	6.44kg	7.75kg	7.89kg
Max Load 200kg	Max Load 200kg	Max Load 200kg	Max Load 200kg
Elevation Angles: 5°, 10°, 15°, 20°	Elevation Angles: 7.5°, 12.5°, 17.5°, 22.5°	Elevation Angles: 5°, 10°, 15°, 20°	Elevation Angles: 7.5°, 12.5°, 17.5°, 22.5°

ATHENA™ SUPINE BREAST SYSTEM

Unmatched adaptability with precision and strength;
that's the power of Athena™.

Included:

Baseplate with Adjustable
Bum Stop
Wingboard with T-bar
and U-grip Handles
Head Support
Mask Adaptor
Arm Cups (CF only)
Upper Hand Grips
Lower Hand Grips
Leg Bolster

Compatible Accessories:

S-type Masks
Breast Masks
Silverman Head Supports
Klarity Cushions™
Vacuum Bags

Powerful
Versatility

Comprehensive hand
placement options
are available with our
U-grip Handles

Arm Cups and
Hand Grips included

Both the Athena
baseplate and the
Athena wingboard are
compatible with all
S-type head masks

Multiple elevation
configurations available


Lower Hand Grips
enable asymmetrical
arm positioning and
offer ergonomic
shoulder suppression

The Athena bum stop is
easily adjustable and locks
securely to prevent movement
and ensure stability

ATHENA™ SUPINE BREAST SYSTEM


OPTIONAL ACCESSORIES	PRODUCT ID
	R322-S424C1-L Klarity White™ Head & Thorax Mask Left Side 3.2mm
	R322-S424C1-R Klarity White™ Head & Thorax Mask Right Side 3.2mm
	R32-2436 White Athena Breast Mask 3.2mm

OPTIONAL ACCESSORIES	PRODUCT ID
	R7661-20C Vacuum Bag for Upper Body Support 58.5x90 cm
	R7626-20C 53cm x 104cm
	R7504-30C 100cm x 70cm
	MR7677-25N Athena Vacuum Bag with mask opening




“We are highly satisfied with the quality and support we get from Klarity. We use Klarity devices across our clinics to ensure the best possible care for our patients.”

— Ainsley Niemkiewicz Fillman, DABR
Director of Physics
Advocate Radiation Oncology
Fort Myers, FL



[KLARITYMEDICAL.COM](https://klaritymedical.com)