

cAutoVerify

Elevate
your workflow



Inspiring excellence in cancer treatment



A more efficient workflow saves time for what really matters

Using an automated workflow to minimize repetitive time-consuming tasks doesn't just save you and your team time, it also increases patient safety. We offer a complete automation solution for both the imaging suite and treatment room, automating patient identity validation and accessory tracking. To streamline the daily workflow within radiation therapy even further, cAutoVerify seamlessly integrates with the Record and Verify system.

Accessory Tracking

Device setup made easy

ALWAYS THE RIGHT ACCESSORIES

With cAutoVerify, you always validate the correct accessories, both patient specific and global.

SEAMLESS INTEGRATION

cAutoVerify automatically validates accessories, reporting directly to the Record and Verify system without any user intervention.

FAST AND RELIABLE

RFID tracking technology brings both speed and reliability to cAutoVerify, and enables completely contactless validation.

ADAPTS TO YOU – NOT THE OTHER WAY AROUND

cAutoVerify adapts seamlessly to your current workflow, streamlining operations, saving time and increasing efficiency whenever the opportunity arises.



Fast device setup



Instant ID validation



Patient Validation

Contactless patient identity validation

ALWAYS THE RIGHT PATIENT

cAutoVerify enhances safety via automation, so you always have the right patient for every treatment.

SEAMLESS INTEGRATION


cAutoVerify automatically validates patient identity, reporting directly to the Record and Verify system without any user intervention.

NEVER FORGETS A FACE

Biometrical facial recognition validates patient identity quickly, without contact, and to the highest level of safety.



Say hello to your **new dream team** in advanced radiation therapy



OUR MISSION IS GLOBAL IN SCOPE

C-RAD's cutting-edge solutions ensure exceptionally high precision, safety and efficiency in advanced radiation therapy, helping to cure more cancer patients and improve their quality of life.

C-RAD – Inspiring excellence in cancer treatment.

C-RAD AB (publ)

C-RAD Positioning AB

Sjukhusvägen 12 K, SE-753 09 Uppsala, Sweden

Telephone +46 18-66 69 30

www.c-rad.com

C-RAD Inc.

70 SE 4th Ave, Delray Beach, FL 33483, USA

Telephone +1 261 742 9260

www.c-rad.com

C-RAD GmbH

Wittestr. 30 K, 13509 Berlin, Germany

Telephone: +49 30 60984 7560

www.c-rad.com

C-RAD (Shanghai)

Medical device Co. Ltd.

Block 2 Room 503,

No 608 Shengxia Road

Pudong new Area, Shanghai, China

